



Repose Flex



Repose<sup>®</sup> Flex

Designed to reduce  
the risk of DVT<sup>1</sup>

repose<sup>®</sup>



Utilising immersion and envelopment, Repose Flex is designed to place the knee in slight flexion (5° to 10°) and allow complete heel offloading while effectively re-distributing the weight of the patient's leg across its surface<sup>2</sup>.

By placing the knee in slight flexion, the risk of Popliteal Vein Compression at the rear of the knee is reduced<sup>1</sup>.

Popliteal Vein Compression is a significant source of Deep Vein Thrombosis for patients undergoing surgery or prolonged bed rest<sup>3</sup>.

Suitable for:

- Patients identified as being at risk of developing a DVT
- Bedbound patients, or patients undergoing surgery
- Protecting patients at all levels of risk of developing all stages of pressure ulcer
- Treatment of pressure ulcers when used as part of a package of care
- Deflatable for convenient storage and transport
- Repose smart valve technology for optimum inflation levels
- Two years warranty as standard



All levels of risk



Easy to clean



Multi-Patient Use

| Description               | Size (L/W/H)               | Ref No. |
|---------------------------|----------------------------|---------|
| Repose Flex - Extra Small | L 420mm x W 700mm x H 90mm | 48030   |
| Repose Flex - Small       | L 460mm x W 700mm x H 90mm | 48031   |
| Repose Flex - Standard    | L 500mm x W 700mm x H 90mm | 48032   |

For more information:

**Call 1300 942 552**

**Email: [info@macmedhealthcare.com](mailto:info@macmedhealthcare.com)**

**Website: [www.macmedhealthcare.com](http://www.macmedhealthcare.com)**

